



As a precaution, gloves and eye protection should be worn when handling drug substances and using test units. Do not puncture any part of unit. Do not remove probe once a capsule is broken. In the unlikely event of a spill, rinse affected area thoroughly with water. Keep out of reach of children. All results are presumptive. If necessary, further testing by a competent drug testing lab can provide additional identification.

Test Preparation	1	Remove the Smart-Tip™ probe from test unit.	Sample suspected substance.	
			Plant Matter	White Powder
			Touch Smart-Tip™ probe to sample. Make sure the tip is completely covered.	Touch Smart-Tip™ probe to sample.

Test Preparation	3	Lock Smart-Tip™ probe into test unit. Press firmly.	STOP Do not continue unless probe is locked below the Flex Pins and cannot be removed from the test unit.	Identification	4	With thumb and forefinger, squeeze chamber 1 until capsule breaks.	5	Gently shake test unit for 2-3 seconds.

Results	The appearance of a yellow color indicates the presence of a Synthetic Cannabinoid , such as:	
	Naphthoylindoles: JWH - 007, 011, 015, 016, 018 , 019, 020, 022 , 072, 073 , 081, 098, 122 , 133 , 180, 182, 200, 210 , 213, 398, 424.	
	Benzoylindoles: AM - 223, 694 , RCS - 4.	
	Nathoylpyrroles: JWH - 030, 031, 145, 147, 307, 309, 368, 369, 370.	
	If no yellow color develops, the above substances are not present. * The test will not react to the following substances:	
	Classical cannabinoids: THC, other constituents of cannabis and their structurally related synthetic analogues, e.g. HU-210, AM-906, AM-411, O-1184.	
	Nonclassical cannabinoids: cyclohexylphenols or 3-arylcyclohexanols such as CP-47, 497-C8, CP-55,940, CP-55,244.	
	Phenylacetylindoles: such as JWH - 203, 250.	
	Proceed to Dilution (steps 6-9), regardless of whether or not a color appeared.	

Dilution	6	With thumb and forefinger, squeeze chamber 2 until capsule breaks.	7	Gently shake test unit for 2-3 seconds.
				Wait 5 minutes. The test window will warm up. This is normal.

Dilution	8	With thumb and forefinger, squeeze chamber 3 until capsule breaks.	9	Gently shake test unit for 2-3 seconds.
				The test window will warm up. This is normal.



Disposal

After Dilution, the test kit is safe to dispose of in the trash.

01-19