



TESTING KITS FOR DRUGS



Procedure Card for IDT3800


KHAT IDENTIFICATION



Substance Test

Read 
Before Use

As a precaution, gloves and eye protection should be worn when handling drug substances and using test units. Do not puncture any part of unit. Do not remove probe once a capsule is broken. In the unlikely event of a spill, rinse affected area thoroughly with water. Keep out of reach of children. All results are presumptive. If necessary, further testing by a competent drug testing lab can provide additional identification.

Test Preparation	1	Remove the Smart-Tip™ probe from test unit.	2	Sample suspected substance.
			Plant Matter Break a small amount of the stem, about the size of the hole (5 mm), and insert into the opening of the kit. 	When possible use a fresh green piece of the stem. Dried pieces may not react as well. Leaves have a weaker reaction.

Test Preparation	3	Lock Smart-Tip™ probe into test unit. Press firmly.	Identification	4	5
		STOP Do not continue unless probe is locked below the Flex Pins and cannot be removed from the test unit.		With thumb and forefinger, squeeze chamber 1 until capsule breaks.	Gently shake test unit for 5 seconds.

Results		The appearance of a pink or red color indicates the presence of Khat .
		If no pink or red color develops, Khat is not present.



Disposal

The test kit is safe to dispose of in the trash.

v2
07/18



Manufactured by IDenta Ltd

www.identa-corp.com

Patented Technology

Made in Israel