



TESTING KITS FOR DRUGS

Procedure Card for IDT5000

AMPHETAMINE/ METHAMPHETAMINE IDENTIFICATION

Substance Test

Read Before Use

As a precaution, gloves and eye protection should be worn when handling drug substances and using test units. Do not puncture any part of unit. Once a capsule is broken, do not resqueeze the chamber or remove the probe. In the unlikely event of a spill, rinse affected area thoroughly with water. Keep out of reach of children. All results are presumptive. If necessary, further testing by a competent drug testing lab can provide additional identification.

Test Preparation	1	Remove the Smart-Tip™ probe from the test unit.	2 Sample the suspected substance.	
			Powder Touch the Smart-Tip™ probe to the sample.	Tablet Firmly press the Smart-Tip™ probe on the tablet while scrubbing in a circular motion. Make sure to completely cover the tip with the substance.

Test Preparation	3	Lock the Smart-Tip™ probe in the test unit. Press firmly.		STOP Do not continue unless the probe is locked in place below the Flex Pins.
		Flex Pins		

Identification	4	With thumb and forefinger , squeeze chamber 1 until the capsule breaks.	5	Gently shake the test unit for 2-3 seconds.

Results		The appearance of a reddish-orange color <u>within one minute</u> indicates the presence of Amphetamine or Methamphetamine .	
		The appearance of a purple color <u>within one minute</u> indicates the presence of an opiate such as Heroin , Morphine or Codeine .	
		The appearance of a dark purple to black color indicates the presence of Ecstasy/MDMA , MDA , MDE or MBDB .	
		If none of the above colors develop, the above listed drugs are not present.	
	Proceed to Dilution (Steps 6-9), regardless of whether or not a color appeared.		

Dilution	6	With thumb and forefinger , squeeze chamber 2 until the capsule breaks.	7	Gently shake the test unit for 2-3 seconds.
				The test window will warm up. This is normal.

Dilution	8	With thumb and forefinger , squeeze chamber 3 until the capsule breaks.	9	Gently shake the test unit for 2-3 seconds.
				The test window will warm up. This is normal.

	Disposal	After Dilution, the test kit is safe to dispose of in the trash.	ENG 01-22
--	-----------------	--	-----------



World Leader in Forensic Products

Manufactured by IDenta Ltd
www.identa-corp.com
Patented Technology
Made in Israel