




Procedure Card for IDT9030



MIX
MULTI DRUGS
IDENTIFICATION
Substance Test



Read Before Use


As a precaution, gloves and eye protection should be worn when handling drug substances and using test units. Do not puncture any part of unit. Once a capsule is broken, do not resqueeze the chamber or remove the probe. In the unlikely event of a spill, rinse affected area thoroughly with water. Keep out of reach of children. All results are presumptive. If necessary, further testing by a competent drug testing lab can provide additional identification.

Test Preparation	1	Remove the Smart-Tip™ probe from the test unit. 	Testing Options	Option 1: If the suspected substance is a powder, proceed to Steps 2-5 on Side 1.
		Option 2 :If the suspected substance is plant material, proceed to Steps 6-11 on Side 2.		

Multi Drug Test

Multi Drugs Test	2	Sample the suspected substance. Powder Touch the Smart-Tip™ probe to the sample. 	3	Lock the Smart-Tip™ probe in the test unit. Press firmly. Flex Pins 	STOP <small>Do not continue unless the probe is locked in place below the Flex Pins.</small>

Multi Drug Identification	4	With thumb and forefinger , squeeze chamber 1 until the capsule breaks. 	5	Gently shake the test unit for 2-3 seconds. 

Multi Drug Results	1	The appearance of a purple-red color indicates that one of these drugs may be present: Heroin, Amphetamine, Methamphetamine, Ecstasy/MDMA, Benzedrone, Buphedrone, Carfentanyl, Cathinone, DMT, EAPB, Fentanyl, Ketamine, mCPP, MDMAI, MDPBP, MDPV, Mephedrone, Methadone, Methcathinone, Methylone, Mexedrone, MPA, PCP, Pentedrone, Pentylone, PMA, PMMA, 2-aminoindane (2-AI), other synthetic/designer drugs.
	2	The appearance of a green or blue-green color indicates the presence of Cocaine or Crack .
	 STOP	If none of the above colors develop, the above listed drugs are not present. Regardless of whether or not any color developed, do not continue to Marijuana testing. STOP and proceed to Disposal.

Disposal The test kit is safe to dispose of in the trash.

ENG 7/13/22



Manufactured by IDenta Ltd
www.identa-corp.com
Patented Technology
Made in Israel






Procedure Card for IDT9030



MIX
MULTI DRUGS
IDENTIFICATION
Substance Test



SIDE 2


Marijuana Test

Marijuana Test	6	Sample the suspected substance.	
		Plant Matter	Hashish
		Touch the Smart-Tip™ probe to the sample. Make sure the tip is completely covered. 	Firmly press the Smart-Tip™ probe on the substance while scrubbing in a circular motion. Make sure to completely cover the tip with the substance. 

Marijuana Test	7	Lock the Smart-Tip™ probe in the test unit. Press firmly. Flex Pins 	STOP <small>Do not continue unless the probe is locked in place below the Flex Pins.</small>

Marijuana Identification	8	With thumb and forefinger , squeeze chamber 2 until the capsule breaks. DO NOT BREAK CHAMBER 1 	9	Gently shake the test unit for 1 minute. 

Marijuana Identification	10	With thumb and forefinger , squeeze chamber 3 until the capsule breaks. 	11	Gently shake the test unit for 10 seconds. 	<small>Wait for the color to fully develop before continuing. This process may take up to 3 minutes.</small>

Marijuana Results	1	The appearance of a blue-violet or purple color indicates the presence of Marijuana/Hashish .
	 STOP	If no blue-violet or purple color develops, Marijuana/Hashish is not present. Regardless of whether or not any color developed, do not continue to Multi Drug testing. STOP and proceed to Disposal.

Disposal The test kit is safe to dispose of in the trash.

ENG 7/13/22



Manufactured by IDenta Ltd
www.identa-corp.com
Patented Technology
Made in Israel