

SWIPPER™

TESTING KITS FOR DRUGS

GHB Test Kits

SW-8000T




**Read
Before
Use**

As a precaution, gloves and eye protection should be worn when handling substances and using test units. Do not puncture any part of unit. Do not remove probe after breaking the ampule. In event of a spill, rinse affected area thoroughly with water. Keep out of reach of children. Store the tubes in a dark place, not exposed to light. All results are presumptive. If necessary, further testing by a competent drug testing lab can provide additional identification.

Test Preparation

1

Remove the SWIPPER™ Cap from the test unit.



Ensure that the test tube remains upright in order to avoid the ampule falling out of the tube.

2

Sample the suspected substance or insert the paper into the liquid.


Touch or Dip



Test Preparation

3

Return the SWIPPER™ Cap to the test tube. Press firmly.



STOP

Do not continue unless the SWIPPER™ Cap is properly secured.

Identification


4

With thumb and forefinger, squeeze the test tube until the ampule breaks.



5

Gently shake test up and down to ensure that the liquid touches the Cap.



Tilt the tube sideways and examine the paper to see if a color change occurs.

Results



Powder

The appearance of a dark blue color indicates the presence of **GHB** (γ-Hydroxybutyric Acid).



Liquid

The appearance of an intense blue color on the paper indicates the presence of **GHB**.

Notes for liquid testing:

- * The test is more reliable for water, tea, black coffee, beer, vodka.
- * The detection limit is approximately 2.5 mg/ml of GHB.
- * The minimum detection limit for wine is approximately 15mg/ml.
- * Water by itself with no GHB will produce a green color.
- * The kit is not reliable for use with milk or any drink with milk, because milk causes a blue color to appear.
- * The kit is not reliable for use with citrus drinks (orange/grapefruit juice).



If no blue color develops or a color other than blue develops, such as green, orange or red, GHB is not present.



Disposal

The test kit is safe to dispose of in the trash.

ENG
v4
01-22



IDenta
World Leader in
Forensic Products

Manufactured by IDenta Ltd
www.identa-corp.com

Patented Technology • Made in Israel

LOT:

EXP: